



Installation and Usage Instructions

Validated Program for Computation of Validation Summary Statistics for Accuracy and Precision

LBA version 1.0 March 12, 2007

CONTENTS:

Application contains two files. One file has the name "**LBA.exe**" contains the files necessary to run the LBA application. The second file with the name "**LBA PDF Files.zip**" contains the validation related documentation to help modify and improve the application as needed for individual users. For examples of the input and output files, see the pdf documents included in the Examples Folder.

INSTALLATION INSTRUCTIONS:

All files necessary to run the LBA application are found in the **LBA.exe** file. This is a self-extracting zip file that will install the application, the template, and the sample data from the publication into a default directory.

The default directory will be:
"C:\Program Files\LBA"

The program directory can be changed, but the files must be kept together, in the same directory, in order for the application to work properly.

The second file can be copied into a user preferred directory following any file copy procedures- COPY and PASTE, Drag and Drop, etc. The individual validation related files are within the zipped file. These files are provided to show the level of validation done on the LBA application and to help users modify and improve the application for their internal needs. Please see the disclaimer below for the use of the validation data. Validation data provided in these files cannot be used by individual users for commercial purposes.

USAGE INSTRUCTIONS:

To run the application, simply launch the "LBA Precision and Accuracy.xls" file and click the "New Summary" button to begin. Detailed usage instructions are provided under "Application Usage Instructions" button found on this page.

Disclaimer:

This software was developed and validated at Procter & Gamble Pharmaceuticals and is being made available free of charge to the scientific community through AAPS. It is not intended for sale or use for any commercial purposes. Neither Procter & Gamble Pharmaceuticals nor AAPS assume any liability for the use of the data generated using this program. It is the responsibility of the user to ensure that the software is appropriate for his/her application and to verify or validate the integrity of the data. Validation data provided in these files cannot be used by individual users for commercial purposes.